

September 26, 2016

Philips Respironics

1001 Murry Ridge Lane
Murrysville, PA, 15668

Re: Airline Travel with the DreamStation Series

To whom it may concern:

The Philips Respironics DreamStation series devices and configurations identified in Table 1 below are in compliance with commercial airline EMI/RFI requirements.

Philips Respironics has designed and tested the identified devices for compliance with section 21, Category M, RTCA DO-160G EMI/RFI requirements as specified in the Code of Federal Regulations 14 CFR 382 *"Nondiscrimination on the Basis of Disability in Air Travel; Final Rule"* and FAA Advisory circular AC No: 91-21B.

Philips Respironics devices approved for airline use have airline use operating instructions included in the User Manual and the following symbol on the label that is located on the bottom of the device:



Complies with
RTCA/DO-160G
section 21, category M

The DreamStation Series User Manuals specify that the device is suitable for use on airlines when the device is operating from an AC or DC power source (Philips Respironics Battery Pack). However, it is not suitable for airline use with any modem, Wi-Fi accessory or humidifier installed.

Table 1 - Devices and configurations tested

- DC power source (Philips Respironics Battery Pack)
- AC power source
- Without heated humidifier
- Without modem or Wi-Fi accessory installed

Part Number	Description
BLX200S15	DreamStation CPAP, BL
BLX200H15	DreamStation CPAP with Humidifier BL
BLX400S15	DreamStation CPAP Pro BL
BLX400H15	DreamStation CPAP Pro with Humidifier BL
BLX500S15	DreamStation Auto CPAP BL
BLX500H15	DreamStation Auto CPAP with Humidifier BL
BLX600S15	DreamStation BiPAP Pro BL
BLX600H15	DreamStation BiPAP Pro with Humidifier BL
BLX700S15	DreamStation Auto BiPAP BL
BLX700H15	DreamStation Auto BiPAP with Humidifier BL
DEX400S13	DreamStation CPAP Pro GER
DEX500S13	DreamStation Auto CPAP GER
DEX700S13	DreamStation Auto BiPAP GER
EEX200S15	DreamStation CPAP EE
EEX200H15	DreamStation CPAP with Humidifier EE
EEX400S15	DreamStation CPAP Pro EE
EEX400H15	DreamStation CPAP Pro with Humidifier EE
EEX500S15	DreamStation Auto CPAP EE
EEX500H15	DreamStation Auto CPAP with Humidifier EE
EEX600S15	DreamStation BiPAP Pro EE
EEX600H15	DreamStation BiPAP Pro with Humidifier EE
EEX700S15	DreamStation Auto BiPAP EE
EEX700H15	DreamStation Auto BiPAP with Humidifier EE
ESX200S15	DreamStation CPAP ES
ESX400S15	DreamStation CPAP Pro ES
ESX500S15	DreamStation Auto CPAP ES
ESX600S15	DreamStation BiPAP Pro ES
ESX700S15	DreamStation Auto BiPAP ES
ESX200H15	DreamStation CPAP with Humidifier ES
ESX400H15	DreamStation CPAP Pro with Humidifier ES
ESX500H15	DreamStation Auto CPAP with Humidifier ES
ESX600H15	DreamStation BiPAP Pro with Humidifier ES
ESX700H15	DreamStation Auto BiPAP with Humidifier ES
EUX200S15	DreamStation CPAP EU
EUX200H15	DreamStation CPAP with Humidifier EU
EUX400S15	DreamStation CPAP Pro EU
EUX400H15	DreamStation CPAP Pro with Humidifier EU
EUX500S15	DreamStation Auto CPAP EU
EUX500H15	DreamStation Auto CPAP with Humidifier EU
EUX600S15	DreamStation BiPAP Pro EE
EUX600H15	DreamStation BiPAP Pro with Humidifier EE

EUX700S15	DreamStation Auto BiPAP EE
EUX700H15	DreamStation Auto BiPAP with Humidifier EE
FRX400S14	DreamStation CPAP Pro w/ C-Flex+ FR
FRX401S14	DreamStation CPAP Pro w/ C-Flex FR
FRX500S14	DreamStation Auto CPAP w/ A-Flex FR
FRX501S14	DreamStation Auto CPAP w/ P-Flex FR
FRX700S14	DreamStation Auto BiPAP with BiFlex FR

GBX200H15	DreamStation CPAP w/Humid, GB
GBX200S15	DreamStation CPAP, GB
GBX400H15	DreamStation CPAP Pro w/Humid, GB
GBX400S15	DreamStation CPAP Pro, GB
GBX500H15	DreamStation Auto CPAP w/Humid, GB
GBX500S15	DreamStation Auto CPAP, GB
GBX600H15	DreamStation BiPAP Pro w/Humid, GB
GBX600S15	DreamStation BiPAP Pro, GB
GBX700H15	DreamStation Auto BiPAP w/Humid, GB
GBX700S15	DreamStation Auto BiPAP, GB

INX200S15	DreamStation CPAP, INTL
INX200H15	DreamStation CPAP w/ Humid, INTL
INX200T15	DreamStation CPAP w/ Humid/HT, INTL
INX400S15	DreamStation CPAP Pro, INTL
INX400H15	DreamStation CPAP Pro w/ Humid, INTL
INX400T15	DreamStation CPAP w/ Humid/HT, INTL
INX500S15	DreamStation Auto CPAP, INTL
INX500H15	DreamStation Auto CPAP w/ Humid, INTL
INX500T15	DreamStation CPAP w/ Humid/HT, INTL
INX600S15	DreamStation BiPAP Pro, INTL
INX600H15	DreamStation BiPAP Pro w/ Humid, INTL
INX600T15	DreamStation BiPAP w/ Humid/HT, INTL
INX700S15	DreamStation Auto BiPAP, INTL
INX700H15	DreamStation Auto BiPAP w/ Humid, INTL
INX700T15	DreamStation BiPAP w/ Humid/HT, INTL

ITX200S15	DreamStation CPAP IT
ITX200H15	DreamStation CPAP with Humidifier IT
ITX400S15	DreamStation CPAP Pro IT
ITX400H15	DreamStation CPAP Pro with Humidifier IT
ITX500S15	DreamStation Auto CPAP IT
ITX500H15	DreamStation Auto CPAP with Humidifier IT
ITX600S15	DreamStation BiPAP Pro IT
ITX600H15	DreamStation BiPAP Pro with Humidifier IT
ITX700S15	DreamStation Auto BiPAP IT
ITX700H15	DreamStation Auto BiPAP with Humidifier IT
ITX400H15W	DreamStation CPAP Pro w/Humid/WiFi, IT
ITX500H15W	DreamStation Auto CPAP w/Humid/WiFi, IT

NDX200S15	DreamStation CPAP ND
NDX200H15	DreamStation CPAP with Humidifier ND

NDX400S15	DreamStation CPAP Pro ND
NDX400H15	DreamStation CPAP Pro with Humidifier ND
NDX500S15	DreamStation Auto CPAP with Humidifier ND
NDX500H15	DreamStation Auto CPAP with Humidifier ND
NDX600S15	DreamStation BiPAP Pro ND
NDX600H15	DreamStation BiPAP Pro w/Humid, ND
NDX700S15	DreamStation Auto BiPAP, ND
NDX700H15	DreamStation Auto BiPAP w/Humid, ND
TRX200S15	DreamStation CPAP TR
TRX200H15	DreamStation CPAP with Humidifier TR
TRX400S15	DreamStation CPAP Pro TR
TRX400H15	DreamStation CPAP Pro with Humidifier TR
TRX500S15	DreamStation Auto CPAP TR
TRX500H15	DreamStation Auto CPAP with Humidifier TR
TRX600S15	DreamStation BiPAP Pro TR
TRX600H15	DreamStation BiPAP Pro with Humidifier TR
TRX700S15	DreamStation Auto BiPAP TR
TRX700H15	DreamStation Auto BiPAP with Humidifier TR
ZAX200S15	DreamStation CPAP ZA
ZAX200H15	DreamStation CPAP with Humidifier ZA
ZAX400S15	DreamStation CPAP Pro ZA
ZAX400H15	DreamStation CPAP Pro with Humidifier ZA
ZAX500S15	DreamStation Auto CPAP ZA
ZAX500H15	DreamStation Auto CPAP with Humidifier ZA
ZAX600S15	DreamStation BiPAP Pro ZA
ZAX600H15	DreamStation BiPAP Pro with Humidifier ZA
ZAX700S15	DreamStation Auto BiPAP ZA
ZAX700H15	DreamStation Auto BiPAP with Humidifier ZA

If you have any other questions regarding our products, please feel free to call the Philips Respironics Customer Services department 1-724-387-4000 or +49 8152 93060. You can also use the following address:

Philips Respironics
1001 Murry Ridge Lane
Murrysville, PA 15668

Respironics Deutschland
Gewerbestrasse 17
82211 Herrsching, Germany



Becky Pishney
Regulatory Affairs Engineer, Sleep Therapy